

ISOMARK Fluid Compounds

ISOMARK Grade / Part Number	Description and Key Features	Working Life @ 25°C	Curing Time @ 25°C
F-1 Grey (F-1G)	Fast curing compound, ideal for cold conditions. Specially formulated to give optimum contrast characteristics for comparator microscopy	1 mins.	7 mins.
F-2 Grey (F-2G)	Slower curing version of F1 Grey. Suitable for ambient temperatures above 20°C.	5 mins.	20 mins
F-1 Black (F-1B)	Fast curing compound. Ideal for cold conditions. Specially formulated for high magnification optical microscopy with incident illumination.	1 min.	7 mins.
F-2 Black (F-2B)	Slower curing version of F1 Black. Suitable for ambient temperatures above 20°C	5 mins.	20 mins.
F-3 Black (F-3B)	Extra slow curing version of F1 Black. Suitable for ambient temperatures above 20°C. Ideal for use with the Isomark spraying equipment.	15 mins.	60 mins.
F-3 White (F-3W)	Extra slow curing version of F1 Black. Suitable for ambient temperatures above 20°C. Ideal for use with the Isomark spraying equipment.	15 mins.	60 mins.
CLEAR (Available as F-1, F-2 and F-3, e.g. F-1 Clear)	For use as a suitable background when removing particulate materials.	As either F-1, F-2 or F-3.	As either F-1, F-2 or F-3.
Footprint	Specially formulated for casting footmarks in mud, sand, clay etc.	30 seconds.	2-3 mins.